मुख्यमंत्री नि:शुल्क जाँच योजना

Rajasthan Medical Services Corporation Limited

Gandhi Block, Swasthya Bhawan, Tilak Marg, C-Scheme, Jaipur - 302005

Ph. No. 0141-2223887, Fax: 0141-2228065

E-mail: edepmrmsc-rj@nic.in

CIN:U24232RJ2011SGC035067

Website: www.rmsc.health.rajasthan.gov.in

F-8(219) RMSC/EPM/M-2/2020-21/NIB-444/ 184

dated: 03 11/2020

Order for Rate Contract

Item/Bidder status: Importer Indicative Quantity: 1

S. No.	Item	Detail	
1.	Name of Equipment / Item	HPLC-Mass spectrophotometer with accessories	
2.	Name & Address of Firm	M/s Macflow Engineering Pvt. Ltd. E-59, 1 st Floor, Pratapganj Road, Pandav Nagar, Delhi-1100991Email: macflow@macflow.in	
3.	Make	Waters	
4.	Model	E2695	
5.	Packing Unit	Each	
6.	R/C No. & Validity	(R.C. No 617 Valid upto 30.09.2022)	
7.	Net Rate Per packing Unit (Rs.)	1,03,46,100.00 (Rs. One Crore Three Lac Forty Six Thousand One Hundred Only) HSN Code 90272000	
8.	RMSC Surchage @5%	Rs. 5,17,305.00	
9.	Rates for Supply to RMSC from Suppliers		
	(i) Basic Rate (In Rs.)		10346100.00
	(ii) IGST@ 18 % (In Rs.)		1862298.00
	(iii) CGST @ 0 % (In Rs.)		0.00
	(iv) SGST @ 0 % (In Rs.)		0.00
	(v) Approved Rate for Suppliers From RMSC {(i)+(ii)+(iii)+(iv)} (In Rs.)		12208398.00
10.	Rate for Supply To Consignee From RMSC		
	(i) Basic Rate (In Rs.)		10346100.00
	(ii) RMSC Surcharge @ 5% on (In Rs.)		517305.00
	(iii) CGST @ 9% {On (i)+(ii)} (In Rs.)		977706.45
	(iv) SGST @ 9% {On (i)+(ii)} (In Rs.)		977706.45
	(v) Approved Rate for Consignee From RMSC {(i)+(ii)+(iii)+(iv)} (In Rs.)		12818817.90
11.	Technical specifications of - HPLC-Mass spectrophotometer with accessories is :-		
	HPLC system consisting of Quaternary Gradient pump, Integrated Auto Injector, Online Degasser, Column Heater, PDA/DAD Detector along with software and Accessories. 1. QUATERNARY GRADIENT PUMP: Principle: Low Pressure Mixing - Serial Dual Piston Automatic Gradient Profile Facility with 9 or more gradient curves Quaternary Pressure Gradient pump Operating Flow Range: 0.001–10 mL/min Flow Accuracy: ± 1.0 % Flow Precision: <0.1% RSD Prop. Accuracy: ±0.5% Precision: ≤ 0.5% SD Pressure Range: 350 bar or more Solvent Vacuum Degassing: Built-in , 4-channels Error detection, Leak detection feature should be available. Compressibility compensation should be automatic and continuous. 2. Auto-sampler Sample Capacity: 90 vails or more Injection Volume: 0.1 μL −100 μL Injection Volume Accuracy: ± 1μl		
	Injection Volume Accuracy: ± 1µ1 Injection Volume Precision: <0 1% of RSD		

मुख्यमंत्री निःशुल्क जाँच योजना

Injection Volume Linearity: > 0.999 RSD < 0.5% at 5 – 90 μL Sample Carry Over: Not more than 0.05% from previous injection Temperature Accuracy Sampler: ± 2°C Pressure Range: 350 bar or more 3. Column Heater: Temperature Control Range: Ambient to 65° C Temperature accuracy : $\pm 0.5^{\circ}$ C Colum heater should hold max 2 columns of 30 cm length. Colum switching valve should be supplied along with the instrument for automated method development. Forced air recirculation Peltier-based column heater/coolermodule DETECTORS: Photo Diode Array/ Diode Array Detector 4. Wave length Range 190-800 nm Diode Element minimum 512/1024 or more Wave length Accuracy ±1 nm Wave length Precision/repeatability ± 0.1 nm Noise Wide slit: $\leq \pm 10 \,\mu\text{AU}$ at 254 nm Response Time = 2 s Slit width - 1 nm Drift < 1 mAU/h at 254 and 520 nm with dry cell/water at 1.0 mL/min Temperature Path length 10 mm (or suitable) Cell Volume Approx. 8-10µL Light source: Suitable for the system Lamp Hour >2000 hr lifetime Max. Data Collection Rate: 80 Hz or more including 3D data acquisition. Lamp optimization software for low noise performance. Single Quadrupole LCMS Detector Mass range (m/z) 30 – 1100 with unit mass resolution Detector: Photo multiplier with high energy conversion Dynode Ion source type: Heated Electrospray Ionization Source type: Atmospheric Pressure Ionization Electrospray Ionization flow rate: Min 1.0 mL/min or higher Supported modes: Full Scan, Single Ion Monitoring (SIM) & simultaneous Full Scan/SIM Scan rate: Min 10,000 (Da/s) or higher Scan Speed: Upto 8 Hz or better SIM sensitivity (positive): 10 pg Polarity switching: max 25 ms orbetter Mass accuracy $\leq \pm 0.2$ Da Mass stability: better than 0.1 Da over 24 hours. Digital dynamic range: ≥ 10⁴ Vacuum interlock should be present. Software: Single software should be able to control LC and MS together Noise free Nitrogen generator along with built in compressor should be supplied. Software: 6. Same software should be able to control all modules of HPLC system. To control, acquisition, online display, processing peak point integration and reporting HPLC data, Full 32/64 Bit Architecture software Software must be 21 CFR part 11 compliant fulfilling all effective regulatory requirements. Windows 10/8 environment or suitable Data reports, online help and wizards Data Integrity, Advanced Security, Audit Trials System suitability min 5 Parameter can be checked Template/Wizards saving & Auto run of templates along with simplified automatic common functions. Calibration curves facility Facility for data security, audit trails and electronic signatures etc., should be available for GLP and 21 CFR compliance. Columns (Pore Size 5μ): (i)C-18 Column: 4 Nos. (Two of 250x4.6 mm & Two of 150 x 4.6 mm).



(i)C-18 Column: 4 Nos. (One of 100x2.1 mm & One of 150 x 2.1 mm).

(ii)C-8 Columns: One of 250x4.6 mm Columns (Pore Size sub 2μ): (major use)

(ii)C-8 Columns: One of 100x2.1 mm



मुख्यमंत्री निःशुल्क जाँच योजना

- 8. Accessories:
 - 1.(A) Vials: 1 ml pack of 100 vials & 2 ml pack of 100 vials (One Each)
 - (B)Low insert vials (for low volumes): A pack of 100 vials
 - 2.Filtration Assembly consisting of
 - (A) Sample Filtration Kit One Nos.
 - (B)Membranes [Type: Dual (Aqueous & Organic solvents)]
 - Size 13 mm diameter with Poresize 0.45 μ 10 Pkt of 100 membranes
 - Size 47 mm of 0.45 µ Poresize; Qty-10 pkt of 100 membranes
 - (C)Pre-filters 10 pkt of 100 circles
 - (D)Solvent filtration kit -One No
 - (E)Imported Oil Free Vacuum pump One No.
 - (F)Nylon Syringe filters (13 mm; 0.45μ) 2 Box of 100 filters.
 - 3. MS Grade solvents of reputed brand (Merck/Fisher)
 - (A) Methanol 5 x2.5 Lit.
 - (B) Acetonitrile- 5 x2.5 Lit
 - (C) Water- 5 x2.5 Lit.
- Suitable PC & Printer with 10 KVA UPS of 60 minutes back up of reputed brands (PC-HP/Dell; Printer Brother/Canon; UPS- IPS/APS/GE) specification as under:

Processor - Intel i7 Latest generation; RAM - 16 GB; Hard disk - 2 TB; Graphic Card; DVD writer; 21" TFT screen; LAN Port; USB 2.0 Ports (4 Nos.); Wi-Fi; Multimedia Keyboard; Optical Mouse

- 10. Operating System Preloaded Windows 10/8
- 11. Printer- Laser Colour Printer with duplex printing and LAN port.
- 12. Documents /Training and other notes:
 - (i) The equipment should be certified by European CE marking or UL marking.
 - (ii)Should have safety certificate from a competent authority like CE/Equivalent or valid detail of electrical and functional safety test. Copy of the certificate/test report shall be produced along with the technical bid
 - (iii)Comprehensive guarantee period with spare parts for at least three years from the date of installation.
 - (iv)Instruction Manuals (Hard Copy or Soft copy as applicable) along with standard Operation Procedures (SOP) and appropriate Dust cover for the instrument.
 - (v)Installation qualification (IQ), Operational Qualification (OQ) and Performance Qualification (PQ) documents should be provided at the time of supply of the instrument.
 - (vi)Satisfactory technical and application training to the personals at the site after installation.
 - (vii)CMC will be given @3% of net rate-inclusive of GST (as applicable) & yearly escalation of 5% on last year's CMC price. The CMC may be awarded for five years (on yearly basis) after completion guarantee period.
 - (viii)Minimum of two routine maintenance services/inspections to be provided annually and any number of breakdown emergency calls should be attended during the guarantee period.
 - (ix)Price quoted should be inclusive of their complete installation in all respects at site as per purchase order.
 - (x)Final technical approval only after demonstration of the quoted product.
 - (xi)Any other perquisite required (if any) with specifications should be included with the instrument to run and complete the installation.

Executive Director (EPM)
Rajasthan Medical Services Corporation Ltd.,

dated: 03 11 2020

F-8(219) RMSC/EPM/M-2/2020-21/NIB-444/ 184

Copy forwarded to the following for information & necessary action:

- 1. P.S. to Hon'ble Health Minister, Medical & Health Deptt. Govt. of Rajasthan.
- 2. P.S. to PHS, Medical Health & FW Deptt. Rajasthan, Jaipur.
- 3. PS to Pr. Secy/Secy, Medical Education Department, Rajasthan, Jaipur
- 4. P.S. to Mission Director NHM Jaipur.
- 5. PA to Managing Director, RMSCL
- 6. Director, (PH)/(Aids)/(FW)/(RCH)/(IEC)/(ESI).
- 7. Principal, Medical College, Ajmer, Bikaner, Jaipur, Jodhpur, Kota, Jhalawar& Udaipur.
- 8. Principal R.U.H.S. College of Dental Science Jaipur.
- 9. Principal, RUHS College of Medical Scinces, Jaipur.
- 10. Superintendent Japuriya Hospital Jaipur (Associated with RUHS College of Medical Scinces, Jaipur).
- 11. Superintendent, SMS Hospital, Zenana Hospital, MahilaChikitsalya, Mental Hospital, J.K. Lon Hospital, T.B & Chest Hospital, Jaipur.
- 12. Superintendent Associated Group of Hospital-Ajmer, M.G. Hospital, Jodhpur, Ummaid Hospital- Jodhpur, PBM Hospital-Bikaner, MBS Hospital-Kota & A.G. Hospital-Udaipur.



मुख्यमंत्री निःशुल्क जाँच योजना

- 13. Director, Medical & Health Services (ESI), Mobile Surgical Unit, Raj. Jaipur.
- 14. Addl. Director, Medical & Health Services (PW)/ (RH)/(Hospital Admn.)
- 15. Joint Director, Medical & Health Services, (Zone) Jaipur, Ajmer, Udaipur, Jodhpur, Bikaner, Kota &Bharatpur.
- 16. All P.M.O.'s.
- 17. All CM&HO's.
- 18. All Dy. CM&HO's.
- 19. AGM (IT), RMSC, Jaipur.
- 20. M/s Macflow Engineering Pvt. Ltd., E-59, 1st Floor, Pratapganj Road, Pandav Nagar, Delhi-1100991Email: macflow@macflow.in

21. Guard file.

Executive Director (EPM)

nan Medical Services Corporation L.

Rajasthan Medical Services Corporation Ltd., Jaipur