



D- Block, Swasthya Bhawan, Tilak Marg, C-Scheme JAIPUR (Raj)

Ph. No. 0141-2223887, Fax No. 0141-2228065

E-Mail –<u>edepmrmsc-ri@nic.in</u>

No. F-8() RMSC/EPM/M-4/Sec. B/MVA Kit/ NIB-53/ 13-14/ 40 70

Dated: 30/7/1

CORRIGENDUM/ADDENDUM

Subject: - Amended date extension and technical Specifications of bid document for the MVA KIT (Syringe with accessories) under NIB No. F-8() RMSC/EPM/M-4/ Sec. B/MVA Kit/ NIB-53/ 13-14/3746 Dated:19.06.2014.

1. In reference to subject cited above NIB No. F-8() RMSC/EPM/M-4/Sec. B/MVA Kit/ NIB-53/ 13-14/3746 Dated: 19.06.2014, the last date for sale of bid, receipt of bid & date of opening of technical bid is hereby extended as below:-

Existing Dates			Extended Dates		
Last date for Sale of Bid Form	Last date of Receipt of Bid Form	Date of Opening of Technical Bid	Last date for sale of Bid Form	Last date of Receipt of Bid Form	Date of Opening of Technical Bid
1.	2.	3.	4.	5.	6.
30-July-2014 5:00 PM	31-July-2014 Upto 1:00 PM	31-July-2014 From 3:00 PM	26-Aug-2014 5:00 PM	27-Aug-2014 Upto 1:00 PM	27-Aug-2014 From 3:00 PM

2. In reference to above cited subject and NIB, the technical specification of items (Annexure-c) are here by amended as below:-

Technical Specification for MVA KIT (Syringe with accessories)

Syringe Type: Double valve manual vacuum aspiration syringe

Material: Plastic Capacity: 50-60 cc.

Pressure: 24-26 inches of Mercury

Plunger: Double arm plunger with "O" Ring Collars stop along with retaining clips.

Valve: Made of Silicon double pinch valve with a liner.

Cannula: Maliable plastic, 4, 5, 6, 7, 8, 9, 10 and 12 mm (total = 8) in size, in 5 different colours.

Compatible with the syringe without an adapter.

Additional three (3) "O" rings & Silicon oil for lubrication, 3 (three) extra collar stop to be supplied with each kit.



Kit to be supplied in transparent plastic box which is also autoclavable.

Sterilization: All parts of the syringe should be cold sterilizable with Gluteraldehyde & 0.5% Chlorine, should be also autoclavable at 121° C.

Reusability: The kit should be reusable at least for 25 times under defined sterilization conditions. A certificate in this regard should be furnished from a Government approved lab.

Shelf life: Minimum 2 yrs. from the date of delivery.

Manufacturer should be ISO 13485 certified.

This is issued with the approval of MD, RMSCL, Jaipur.

Executive Director (EPM)

RMSC, Jaipur